

INSTRUCTION MANUAL

Manifold and Filtration System

Model No.

MultiVac 310 / 610 - MS

MultiVac 300 / 600 - MS

MultiVac 300 / 301 / 600 / 601 - MB

MultiVac 310 / 610 - MS - T / A

MultiVac 301 / 601 - MB - T / A

Tech support



ROCKER

Please read this instruction manual before using this product.

Table Of Contents

1. Important Notice	1
2. Unpacking	2
(1) Manifold	2
(2) Vacuum Filtration System	4
3. Main Part Diagram	6
(1) Manifold	6
(2) Vacuum Pump	6
(3) Parts Material List	7
4. Installation and Operation	8
(1) Preparation and Installation of Manifold	8
(2) Vacuum Filtration System	9
2.1 System with Alligator 200 Diaphragm Liquid Pump	9
2.2 System with Rocker 400 Vacuum Pump	10
5. Maintenance	11
6. Troubleshooting	12

1. Important Notice


This instrument is designed for laboratory usage only. Please read this manual carefully before installing and operating. The instrument shall not be modified in any way. Any modification will void the warranty and may result in potential hazard. We are not responsible for any injury or damage caused by any non-intended purposes and modifying the instrument without authorization.

1. Check the voltage specified on the name plate and ensure it matches the line voltage in your location.
2. Install the instrument in a clean, dust-less and ventilated area under 40°C.
3. Never use the pump with any chemical, corrosive, flammable or toxic material.
4. Before operation, please check the compatibility of the contact materials and used medium.
5. Please use liquid pump with manifold and membrane to prevent clogging.
6. The protection device of Rocker 400 vacuum pump:
 - To ensure proper function of the Prote™ protection device, the pump must be installing on a flat level surface and should not be moved when operation. If the buoy closes during operation, it must be manually reset by uninstalling and reinstalling the filter cartridge.
 - The filter cartridge is used to absorb moisture and dust. Please replace the cartridge to maintain pumping efficiency when it becomes dirty or color-changing.
 - The pump has a thermal protection device that automatically shuts off when overheated.
7. Do not use any lubricant, which may damage the pump.
8. In case of any damage, please contact the manufacturer or your service agent for service to avoid hazard.
9. Please discard packing material according to local related regulations.
10. Operating condition
 - (a) Ambient temperature: 5~ 40 °C
 - (b) Relative humidity: 80% RH Max.
 - (c) Power supply: (Alligator 200) 100-240V~, 50/60Hz
(Rocker 400) 100-120V~, 50/60Hz or 200-240V~, 50/60Hz
 - (d) Altitude: up to 2000 m
 - (e) Pollution degree: II
 - (f) Indoor use

2. Unpacking

Please check if the package is complete without any damage before unpacking. When unpacking, please make sure you have all accessories that indicated on the list. If there is any problem, please keep the serial number along with packing case and contact your local distributor immediately for assistance.

(1). Manifold



MultiVac 300 / 310 - MS



MultiVac 600 / 610 - MS



MultiVac 300 - MB



MultiVac 600 - MB

	Standard Package Includes:		
	Branches / Material	SF Adaptor	Stopper Adaptor
MultiVac 310 - MS	3 Places Aluminum Manifold		V
MultiVac 610 - MS	6 Places Aluminum Manifold		V
MultiVac 300 - MS	3 Places Stainless Steel Manifold		V
MultiVac 600 - MS	6 Places Stainless Steel Manifold		V
MultiVac 300 - MB	3 Places Stainless Steel Manifold	V	
MultiVac 600 - MB	6 Places Stainless Steel Manifold	V	

• An allen key is included with each manifold..



MultiVac 301 - MB



MultiVac 601 - MB

Model	Standard Package Includes:		
	Branches / Material	SF Adaptor	Stainless Steel Funnel
MultiVac 301 - MB	3 Places Stainless Steel Manifold	V	100 mL x 3
MultiVac 601 - MB	6 Places Stainless Steel Manifold	V	100 mL x 6

- An allen key is included with each manifold.

(2). Vacuum Filtration System

• Manifold with Alligator 200 Diaphragm Liquid Pump



MultiVac 301 / 601 - MB - A

Model	Standard Package Includes:	
MultiVac 301 - MB - A MultiVac 601 - MB - A	MultiVac 301 / 601 - MB, 3 / 6 Places Stainless Steel Manifold	x 1
	Alligator 200, Diaphragm Liquid Pump	x 1
	Dragon 100, Torch Burner	x 1
	MCE Membrane, 47 mm, 0,45 µm	200 pcs
	Silicone Tube	x 2
	Instruction Manual	x 1

- An allen key is included with each manifold.



MultiVac 310 / 610 - MS - A

Model	Standard Package Includes:	
MultiVac 310 - MS - A MultiVac 610 - MS - A	MultiVac 310 / 610 - MS, 3 / 6 Places Aluminum Manifold	x 1
	MF 3, Magnetic Filter Holder (300 mL)	x 3 / 6
	Alligator 200, Diaphragm Liquid Pump	x 1
	Silicone Tube	x 2
	Instruction Manual	x 1

- An allen key is included with each manifold.

• **Manifold with Rocker 400 Vacuum Pump**



MultiVac 301 / 601 - MB - T

Model	Standard Package Includes:	
MultiVac 301 - MB - T MultiVac 601 - MB - T	MultiVac 301 / 601 - MB, 3 / 6 Places Stainless Steel Manifold	x 1
	Rocker 400, Oil Free Vacuum Pump	x 1
	4000 mL PP Vacuum Bottle with drain	x 1
	Dragon 100, Torch Burner	x 1
	MCE Membrane, 47 mm, 0.45 µm	200 pcs
	Silicone Tube	x 2
	Instruction Manual	x 1

• An allen key is included with each manifold.



MultiVac 310 / 610 - MS - T


Model	Standard Package Includes:	
MultiVac 310 - MS - T MultiVac 610 - MS - T	MultiVac 310 / 610 - MS, 3 / 6 Places Aluminum Manifold	x 1
	MF 3, Magnetic Filter Holder (300 mL)	x 3 / 6
	Rocker 400, Oil Free Vacuum Pump	x 1
	4000 mL PP Vacuum Bottle with drain	x 1
	Silicone Tube	x 2
	Instruction Manual	x 1

• An allen key is included with each manifold.


3. Main Part Diagram

(1). Manifold

Stainless Steel Manifold



Aluminum Manifold



Position	Designation	Position	Designation
1	Control Valve (Knob)	4	Quick-Fit Hose Barb
2	Handle	5	Funnel Adaptor
3	Adjustable Vent Valve		

- Adjustable vent valve exits only on stainless steel manifold.

(2). Vacuum Pump

Rocker 400 Oil Free Vacuum Pump



Alligator 200 Diaphragm Liquid Pump



Position	Designation	Position	Designation
1	Power Switch	6	Vacuum Gauge
2	Air Outlet	7	Moisture Trap
3	Air Inlet	8	Vacuum Regulator
4	Liquid Inlet (IN)	9	Prote™ Filter Cartridge
5	Liquid Outlet (OUT)	10	Power Adaptor

(3). Parts Material List

- MultiVac, Manifold

Models	Material contacts with filtrates
MultiVac 310 / 610	Aluminum, Viton
MultiVac 300 / 600	SS316, Viton

- Alligator 200, Diaphragm Liquid Pump



Part	Material contacts with filtrates
Wetted parts	SS316, FKM, Silicone, PP, PPS

- Rocker 400, Oil Free Vacuum Pump


Part	Material	Part	Material
Pump Head	Aluminum alloys	Piston Seal	PTFE complex
Valve Plate(s)	FKM		

4. Installation and Operation

(1). Preparation and Installation of Manifold

- A. Make sure all control valves are off  before use, rotate the control valves counterclockwise to vertical position  for filtration.
- B. Change the funnel adaptor to fit the desired funnel or filter holder if necessary.
- Various kind of funnel adaptors are only compatible for **Stainless Steel Manifold**.
 - **How to change the funnel adaptor(s) and remove the membrane support**
 1. Loosen the fixing screw on the rear of control valve by allen key (2 mm).
 2. Unscrew the funnel adaptor counterclockwise and remove the membrane support.
 3. Replace it with the desired funnel adaptor and lock it back in a clockwise direction.
 4. Tighten the fixing screw to prevent funnel adaptor from loosening.

Installation of Various Filter Holder and Funnel



1. No. 8 Stopper Funnel Adaptor

- Compatible with various kinds of filter holders with No. 8 silicone stopper, such as plastic and stainless steel filter holder and most disposable filter funnels.

2. ROCKER SF Funnel Adaptor (Stainless Steel Manifold only)

- Fitting to (left) ROCKER SF 100 / 300 / 500 mL funnel and (right) Pall, Sentino® 4870 (100 mL), 4871 (250 mL) disposable plastic funnel.

3. SS316 Adaptor for Millipore Disposable Funnel (Stainless Steel Manifold only)

- Fitting to Millipore, Microfil® disposable plastic funnel 100 / 250 mL.

4. SS316 Adaptor for Sartorius Disposable Funnel (Stainless Steel Manifold only)

- Fitting to Sartorius, Microsart® disposable plastic funnel 16A07.



- C. Connect installed MultiVac manifold to either vacuum pump or liquid pump for filtration.
 - Please ensure to install a vacuum bottle between the manifold and vacuum pump to capture the waste liquid and prevent liquid from being sucked into the vacuum pump.
 - The inlet and outlet direction of MultiVac manifold is interchangeable. Please disassembled inlet and outlet by hands or tool and change to desired direction.
- D. The adjustable vent valve can help drain residual liquid in tube and avoid “Water Hammer”.
 - Turn the vent valve counterclockwise to drain the remaining liquid out if necessary.

(2). Vacuum Filtration System Operation

2.1 Vacuum Filtration System (With Alligator 200 Diaphragm Liquid Pump)



Filtration System Assembly Diagram

- A. Connect the outlet of installed manifold to inlet (IN) of Alligator 200 liquid pump together with silicone tube as shown above.
 - The tube length between manifold and liquid pump should be as short as possible.
 - The length of silicone tube can be cut as needed.
- B. Connect a long tube to outlet (OUT) of Alligator 200 and put the other side of tube into sink or a container to drain and collect the waste liquid.
- C. Install the appropriate membrane in filter holder or funnel.
- D. Activate on the pump and pour the sample liquid into the funnel and open the control valve to “vertical” position  to start filtration.
- E. After finishing filtration, close the control valve  and retrieve the membrane.
- F. SF stainless steel funnel can be fast flame-sterilized by Dragon 100 torch burner or autoclaved before proceeding next filtration.

*** Warning:**

(1) Do NOT use cable tie or any similar items to secure the pump's outlet and silicone tube together to prevent any bursting and potential damage to the pump.



(2) Do NOT bend silicone tube or squeeze it to avoid poor drainage and potential damage.

(3) Do NOT pump dirty or impure solution without filtration to avoid clogging or damaging.

2.2 Vacuum Filtration System (With Rocker 400 Vacuum Pump)







Filtration System Assembly Diagram

- A. Connect the installed manifold, vacuum bottle and Rocker 400 vacuum pump with 2 silicone tubes as shown above.
- B. Install the appropriate membrane in filter holder or funnel.
- C. Switch on the pump and pour the sample liquid into the funnel and open the control valve to “vertical” position  to start filtration.
 - Please monitor the liquid level in waste bottle to prevent liquid from entering the pump.
- D. After finishing filtration, close the control valve  and retrieve the membrane.
- E. ROCKER SF stainless steel funnel can be fast flame-sterilized by Dragon 100 torch burner or autoclaved before proceeding next filtration.

***Warning: Please ensure the liquid level of moisture trap and vacuum bottle do not exceed safety level. Failure to comply could result in serious damage to the pump and void the warranty. If it occurs, please contact the distributor immediately for service.**

5. Maintenance

1. The pump is not autoclavable. Please clean the surface by pure water or 75% ethanol.
2. While doing filtration with Alligator 200 liquid pump, please pump the pure water to clean the tubing after every use to prevent any dirty is left over. Keep pumping the air for at least 2 minutes to eliminate the liquid.
3. While doing filtration with Rocker vacuum pump, please remove the liquid in the moisture trap and keep pumping the air for at least 2 minutes to withdraw the residual steam.
 - Please refer to manual of Rocker Series Vacuum Pump to see how to reset the Prote™ protection device and replace cartridge.
4. MultiVac manifold, filter holder and vacuum bottle are autoclavable (121°C, 1 bar, 20 min). Please make sure to rinse manifold with pure water and keep all control valves and vent valve at "OPEN" position to make a complete sterilization. You can also clean the surface of manifold with pure water or 75% ethanol.
5. Funnels, support base and membrane support of stainless steel manifold are made of SS316 which can be fast sterilized by flame or autoclaved.
6. Valve plate(s), O-ring(s), silicone tube(s), muffler and diaphragm are consumables, it is recommended to replace it on a yearly basis or as needed to ensure good operation.

O-ring(s)			
Funnel Adaptor	Control Valve	Plug / Vent Valve	Outlet
			

6. Troubleshooting

Problem	Reason and Solution
Failed to start	<ul style="list-style-type: none"> • Loose plug → Reconnect plug to socket firmly. • Wrong voltage → Reconnect to power specified on name plate. • Vacuum exists → Release vacuum and restart the pump. • Switch is broken → Contact distributor for assistance.
Poor vacuum	<ul style="list-style-type: none"> • O-ring of funnel adaptor or vent valve is broken → Replace the O-ring. • Open vent valve → Close the vent valve. • Outlet is blocked → Clean or stretch the outlet tube. • Improperly vacuum setting → Adjust the vacuum regulator. • Clogged muffler of Rocker 400 pump → Replace a new muffler. • Air leakage → Replace or tighten the tubing and accessories that possible leaks. • Internal wiring loose → Contact distributor for assistance.
Low filtration speed	<ul style="list-style-type: none"> • O-ring of funnel adaptor or vent valve is broken → Replace the O-ring • Open vent valve → Close the vent valve. • Air leakage → Replace or tighten the tubing and accessories that possible leaks. • Unsuitable membrane → Change an appropriate membrane. • Membrane is blocked → Replace a new membrane. • Improperly vacuum setting → Adjust the vacuum regulator. • Silicone tube is bent too much → Replace or stretch silicone tube. • Internal wiring loose → Contact distributor for assistance.

Ordering Information

180301-02	MultiVac 301-MB, 3 Places SS Manifold with ROCKER SF funnel adaptor, 100 mL SF funnel x 3
180601-02	MultiVac 601-MB, 6 Places SS Manifold with ROCKER SF funnel adaptor, 100 mL SF funnel x 6
180300-01	MultiVac 300-MS, 3 Places SS Manifold with No.8 stopper funnel adaptor
180600-01	MultiVac 600-MS, 6 Places SS Manifold with No.8 stopper funnel adaptor
180310-01	MultiVac 310-MS, 3 Places Aluminum Manifold with No.8 stopper funnel adaptor
180610-01	MultiVac 610-MS, 6 Places Aluminum Manifold with No.8 stopper funnel adaptor
190100-71	SS316 Silicone stopper funnel adaptor (long)
190200-20	SS316 Adaptor for ROCKER SF funnel and Pall 4870 / 4871 plastic funnel
190100-73	SS316 Adaptor for Millipore MIHAWG100 / MIHAWG250 plastic funnel
190100-74	SS316 Adaptor for Sartorius 16A07 plastic funnel
180100-01/03/05	ROCKER SF stainless steel funnel 100 mL / 300 mL / 500 mL
180100-31	Stainless steel support screen for ROCKER SF funnel adaptor (#190200-20)
167100-31	Stainless steel support screen for Millipore / Sartorius funnel adaptor
181200-01(02)	Alligator 200, Diaphragm Liquid Pump, AC100-240V, 50/60Hz, US Plug (EU plug)
167400-11(22)	Rocker 400, Vacuum Pump, AC110V,60Hz (AC 220V, 50Hz)
167301-07	Prote™ filter cartridge for Rocker 300 / 400 / 410 / 430
167200-48	4000 mL Heavy Duty PP Vacuum Bottle with drain
200300-01(-P)	MF 3 (Pro), Magnetic Filter Holder (300 mL)
200500-01(-P)	MF 5 (Pro), Magnetic Filter Holder (500 mL)
167120-31	300 mL Glass Filter Holder for VF 6 / VF 12
180300-68	Silicone tube, Ø 8 x 14 mm, 200 cm
167200-38	Silicone tube, Ø 6 x 12 mm, 100 cm
177100-00	Dragon 100, Torch Burner
167100-20	PALL A/E glass fiber membrane filter, 47 mm, 1 µm, 100/PK
167100-52-1	Whatman ME 25 (MCE) membrane, 47 mm, 0.45 µm, black grid, sterile, 100/PK

Rocker Scientific Co., Ltd.

Tel: +886-2-26033311 Fax: +886-2-26036622

E-mail: export@rocker.com.tw <https://www.rocker.com.tw>