

Aluminum Multi-Branch Manifold

Features

- ▶ **Quick-assembly and disassembly**
Each component can be easily removed for cleaning and autoclaving .
- ▶ **Light and handy aluminum**
All parts of MultiVac 310 / 610 manifold are made of aluminum which is light, cost-effective and autoclavable.
- ▶ **Quick-Fit hose barb**
Detachable hose barb can be connected to either side of the manifold quickly towards the vacuum source.
- ▶ **Individual control valve**
Each branch permits individual control.
- ▶ **Compatible with different kinds of filter holders**
MultiVac 310 / 610 - MS manifold includes No.8 silicone stoppers which can accept various types of filter holders at option.

Applications

- ▶ Suspended solid test
- ▶ Microbiological test
- ▶ Vacuum filtration

Ordering Information

180310-01

MultiVac 310 - MS, 3 Places Aluminum Manifold with No. 8 stopper funnel adaptor

180610-01

MultiVac 610 - MS, 6 Places Aluminum Manifold with No. 8 stopper funnel adaptor

200300-01

MF 3, 300 ml Magnetic Filter Holder with lid kit

167110-16

No.8 Silicone Stopper



Specification


Model	MultiVac 310 - MS	MultiVac 610 - MS
Number of branch	3	6
Filter diameter	Optional	Optional
Funnel capacity	Optional	Optional
Hose barb	ID8	ID8
Dimension (LxWxH)	48 x 13 x 16 cm	85 x 13 x 16 cm
Part	Material	Material
Body	Aluminum	Aluminum
Control valve	Aluminum	Aluminum
Adaptor	Aluminum	Aluminum
Hose barb	Aluminum	Aluminum
Handle	Aluminum	Aluminum

How to order?

You can select Combination A, B or C

A

Step 1. Choose the manifold



MultiVac 300 - MB
3 Places Stainless Steel Manifold
(180300-02)



MultiVac 600 - MB
6 Places Stainless Steel Manifold
(180600-02)



Step 2. Choose the SF funnel



Rocker, 100 ml
SF Stainless Steel Funnel
(180100-01)



Rocker, 300 ml
SF Stainless Steel Funnel
(180100-03)



Rocker, 500 ml
SF Stainless Steel Funnel
(180100-05)



Stainless Steel Cover
for SF 100 ml funnel
(180100-11)



Stainless Steel Cover
for SF 300 ml / 500 ml funnel
(180100-13)

B

Step 1. Choose the manifold



MultiVac 300 - MS
3 Places Stainless Steel Manifold
(180300-01)



MultiVac 600 - MS
6 Places Stainless Steel Manifold
(180600-01)



Step 2. Choose the filter holder



Rocker, 300 ml
MF 3, Magnetic Filter Holder
(200300-01)



Rocker, 500 ml
MF 5, Magnetic Filter Holder
(200500-01)



Rocker, 300 ml
VF 6, Glass Filter Holder
(167120-31)


to complete your order

C

Step 1. Choose the manifold



MultiVac 300
3 Places Stainless Steel Manifold Base
(180300-00)



MultiVac 600
6 Places Stainless Steel Manifold Base
(180600-00)



Step 2. Choose the disposable funnel and adaptor



Sartorius, 100 ml / 250 ml
Microsart® disposable plastic funnel
(16A07-10-----N/
16A07-25-----N)



Millipore, 100 ml / 250 ml
Microfl® disposable plastic funnel
(MIHAWG100 / MIHAWG250)



Pall, 100 ml / 250 ml
Sentino® disposable plastic funnel
(4870 / 4871)

+



SS316 Adaptor for Sartorius
disposable plastic funnel
(190100-74)

+



SS316 Adaptor for Millipore
disposable plastic funnel
(190100-73)

+



SS316 Adaptor for Pall
disposable plastic funnel
(190100-75)

* SS316 adaptor is only available for above-mentioned disposable funnels.

Customized Design

We can offer customized adaptor for your funnel

