

BioDolphin Suction Kit



Features

- ▶ **Quick Connection**
Quick connection between handle and adaptor can make any fitting adaptors replaced quickly.
- ▶ **Single-handed tip ejection**
The innovative tip ejection design enables tip ejection by single hand easily.
- ▶ **User-friendly lock switch design**
BioDolphin's ergonomic handle and friendly designed lock switch can keep continuous suction without pressing button all the time.
- ▶ **Multi-functional storage rack**
Exclusive designed rack can accommodate all the adaptors in order. A built-in liquid collection tray can receive liquid residue dripping from the adaptors.
- ▶ **Autoclavable**
The complete kit except rack are made of POM or SS304 material available for autoclave.


Applications

- ▶ Cell culture aspiration
- ▶ Laboratory waste suction

Ordering Information

199100-00 BioDolphin Suction Kit

199200-00 Plus Suction Kit



Standard accessory :

- | | |
|--|--|
| 1. 199100-30 Handle | 6. 199100-49 8-channel SS adaptor (40 mm) |
| 2. 199100-46 1-channel tip adaptor (50 mm) | 7. 199100-75 1-channel tip adaptor with ejector |
| 3. 199100-76 1-channel tip adaptor (150 mm) | 8. 199100-50 8-channel tip adaptor with ejector |
| 4. 199100-47 1-channel SS adaptor (40 mm) | 9. 199100-61 Rack |
| 5. 199100-48 1-channel SS adaptor (80 mm) | |

Specification

| Model | BioDolphin | Plus |
|----------------------|---|--|
| Operation | Press to suck | |
| Hose barb | ID6 | |
| Suction kit material | POM, SS304 | POM |
| Rack material | ABS | - |
| Dimension (LxWxH) | 15 x 15.2 x 19.5 cm | 15 x 7.5 x 5.5 cm |
| Suction kit | 1. Handle 2. 1-channel tip adaptor (50 mm) 3. 1-channel tip adaptor (150 mm) 4. 1-channel SS adaptor (40 mm) 5. 1-channel SS adaptor (80 mm) 6. 8-channel SS adaptor (40 mm) 7. 1-channel tip adaptor with ejector 8. 8-channel tip adaptor with ejector 9. Rack | 1. Handle 3. 1-channel tip adaptor (150 mm) 7. 1-channel tip adaptor with ejector 8. 8-channel tip adaptor with ejector |