# **Solvent Purification System**

## **Features**

#### ► No air pollution, maintenance free

Rocker 300C / Chemker 300 / Lafil 300C vacuum pumps are driven by diaphragm, no need of lubricant, regular oil changes and maintenance; with no oil pollution.

#### Good chemical resistance

Rocker 300C / Chemker 300 / Lafil 300C-VF12 systems all can offer good resistance to a wide range of chemical gases except strong acid and alkali.

#### One-step filtration

VF 12 can filtrate and store solution directly in GL45 storage bottle without transferring from the flask and hence preventing cross-contamination.

#### Suitable for all GL45 storage bottle

Unique filtration adaptor is designed to be compatible with any GL45 storage bottle (1000 ml, 2000 ml, 5000 ml, etc.)

#### International certification

CE certification
IP 30 protection (Rocker 300C)
RoHS certification (Rocker 300C)

# One-Year Warranty

(Rocker 300C / Chemker 300 / Lafil 300C)

# **Applications**

- Solvent purification
- ▶ HPLC mobile phase purification
- ► Vacuum filtration

# **Ordering Information**

# 189302-11(22)

Rocker 300C - VF 12, Solvent Purification System, complete set includes:

- \* Rocker 300C, Vacuum Pump
- \* VF 12, Filtration Set
- \* Silicone Tube (1m)

## 169302-11(22)

Chemker 300 - VF 12, Solvent Purification System, complete set includes :

- \* Chemker 300, Vacuum Pump
- \* VF 12, Filtration Set
- \* Silicone Tube (1m)

### 197331-11(22)

Lafil 300C - VF 12, Solvent Purification System, complete set includes:

- \* Lafil 300C, Vacuum Pump
- \* VF 12, Filtration Set
- \* Silicone Tube (1m)









# **Specification**

Model	Rocker 300C - VF 12		Chemker 300 - VF 12		Lafil 300C - VF 12	
Voltage / Freq.	110V/60Hz	220V/50Hz	110V/60Hz	220V/50Hz	110V/60Hz	220V/50Hz
Max. vacuum	145 mbar		105 mbar		145 mbar	
Max. flow rate	23 L/min	20 L/min	20 L/min	18 L/min	23 L/min	20 L/min
Vacuum regulator	Yes		-		Yes	
Hose barb	ID8					
Model	VF 12					
Funnel	300 ml					
Storage bottle	1000 ml					
Filter diameter	47 mm					