


INSTRUCTION MANUAL

A. Introduction

The Suction cap for GL45 bottle can make a typical tissue bottle work as a suction bottle. It is autoclavable, reusable, easily assembled and ideal for general waste liquid suction.

B. Diagram



See Fig.2: Connect "LIQ" (Liquid inlet) to suction kit and connect "AIR" (air outlet) to vacuum source.

C. Order Information

Cat. No.	Description
197000-60	Suction cap for GL45 bottle
197000-60-500	500ml storage bottle with GL45 suction cap
197000-63	Silicone O-ring for GL45 storage bottle
199100-68-5	Silicone tube (Φ4.2xΦ8mm), 5cm
197110-01	Lafil 100 Plus portable suction system
197403-11(22)	Lafil 400 - BioDolphin portable suction system
199100-00	BioDolphin suction kit (light blue color)