

Filtration Accessories



PES / PC Vacuum Bottle

Specification	PES vacuum bottle, 1200 mL	PC vacuum bottle, 1200 mL
Catalog number	197000-11	197000-11-PC
Graduated increment	100 mL	100 mL
Overflow protection	Yes	Yes
Cap with suction fitting	Yes	Yes
Hose barb	ID8	ID8
Autoclavable	Yes	Yes



PP Heavy-Duty Vacuum Bottle

Specification	PP vacuum bottle, 4000 mL	PP vacuum bottle with drain, 4000 mL
Catalog number	197200-47	197200-48
Connector	Hose barb	Hose barb
Graduated increment	500 mL	500 mL
Overflow protection	Yes	Yes
Hose barb	ID4 ~ ID8	ID4 ~ ID8
Drain device	-	Yes
Autoclavable	Yes	Yes



Multilayer Filtration Adaptor

Specification	Multilayer filtration adaptor
Catalog number	167230-05
Main material	PTFE
Material of membrane support	SS screen
Filter diameter	47 mm
Filter holder	(Optional) Glass filter holder 300 mL



Hoop Holder **NEW**

Specification	Hoop Holder	Hoop Holder
Catalog number	167300-37	189300-37
Applicable pump	Rocker 300 / 400 / 410 Rocker 300C / 400C / 410C	Rocker 300 / 400 / 410 Rocker 300C / 400C / 410C
Applicable filtration apparatus	LF 30 / SF 10 / SF 30	VF 12
Applicable apparatus diameter	Max. 118 mm	Max. 103 mm

* PES : Polyethersulfone, PC : Polycarbonate, PP : Polypropylene